

Dissolved Oxygen and Biochemical Oxygen Demand Analyses



Prepared By
Michigan Department of Environmental Quality
Operator Training and Certification Unit

Note: These Procedures are available in the OTCU Laboratory manual which is available on the OTCU website.

Dissolved Oxygen

D.O.

Amount of “FREE” Oxygen

O_2

In the Water

Required By :

FISH

MICROORGANISMS

(Bacteria)

RESPIRATION

Organic Matter Converted to:

Carbon Dioxide
and
Water



RESPIRATION



Reason for Treatment

Basis of

Secondary Treatment

Basis of BOD Analysis


Nitrification

BIOLOGICAL Oxidation of Nitrogen

From to to
AMMONIA NITRITE NITRATE
(NH₃) (NO₂) (NO₃)



NITRIFICATION



May be Required:

Oxygen Demand Reduction

Ammonia Removal (toxic)

Reason for D.O. Analysis

Used in BOD Analysis

Treatment Processes

Treatment Plant Effluents

Receiving Stream

Maximum D.O. Concentrations in Water (Saturation)

Temperature °C	Max. Concentration mg/L
0	14.6
4	13.1
8	11.9
12	10.8
16	10.0
20	9.2
24	8.5
28	7.9
30	7.6

SAMPLE COLLECTION

Sample


0 to 10 mg/L

Atmosphere


21 %

210,000 mg/L


SAMPLE COLLECTION



SAMPLE COLLECTION



Do Not AERATE Sample



**Do Not Splash Sample in Air
or
Let Air Circulate in Sample**


D.O. Procedure

1888

WINKLER


DISSOLVED
OXYGEN


DETERMINATION



WINKLER Dissolved Oxygen **PROCEDURE**
with the
ALSTERBERG
AZIDE
MODIFICATION

AZIDE Destroys NITRITE





Winkler Method

Iodometric Method



1. Takes Mixing

2. Takes Time

3. Free IODINE released in relation to D.O. in Sample

Winkler Procedure




**1 mL
Manganous
Sulfate**

**1 mL
Alkaline
Azide Solution**



Stopper



Mix Well




**Allow
Floc
to Settle**



Repeat Mixing

**Settle
Again**



Contact and Time



White - No Dissolved Oxygen



**1 mL
Sulfuric
Acid**


**Sulfamic
Acid**



Mix

**Measure
100 mL**

**Iodine
Solution**




Titration

**Sodium Thiosulfate
(Thio)**

OR

**Phenylarsene Oxide
(PAO)**



End Point Indicator





Starch



Blue


Clear





Mg/L D.O. =

$$\frac{\text{mL of Titrant} \times \text{Normality of Titrant} \times 8 \times 1000}{\text{mL of Sample}}$$



$$\frac{0.0125 \text{ N} \times 8 \times 1000}{100} = 1$$

Mg/L D.O. =

mL of Titrant X 1

When 100 ml of Sample is Titrated with 0.0125 Normal Titrant

Each mL used Equals 1 mg/L of D.O.

If

4.6 mL

is Used

Then the D.O. in the Sample is

4.6 mg/L

1 mL of Titrant = 1 mg/L D.O.

**With Any of the Following Combinations of
Sample Volumes and Normalities**

Sample Volume

Normality

100 mL

0.0125 N

200 mL

0.0250 N

300 mL

0.0375 N

Outline Of Winkler Dissolved Oxygen Procedure

3.1

Add
1 mL
 MnSO_4 Soln.
and
1 mL
Alkali-iodide-azide
Reagent




Yellow
To
Brown
Floc,
D.O. Present



White
Floc,
No D.O.


3.11

Mix By
Inverting
and
Allow To
Settle



3.12

Repeat
Mixing
and
Settling




3.13

Add
1 mL
 H_2SO_4
and
Mix




Carefully
Collect
Sample
In 300 mL
BOD Bottle




Titration of Iodine Solution

3.21

Pour
100 mL
Into Flask




Reddish-Brown



3.22


Titrate
With
THIO



Pale
Yellow

3.23


Add
Starch
Indicator



Blue

3.24



Titrate
to
Clear



Clear

Winkler Method

Iodometric Method



1. Takes Mixing

2. Takes Time

3. Free IODINE released in relation to D.O. in Sample

D.O. Procedure


Electrode Methods

MEMBRANE ELECTRODE METHOD



The membrane electrode is composed of two solid metal electrodes in contact with supporting electrolyte separated from the test solution by a gas permeable membrane. Oxygen dissolved in the sample diffuses through the membrane on the DO probe and is chemically reduced (accepts electrons), producing an electrical current between the anode and cathode in the probe. The amount of current is proportional to the concentration of DO. Following proper calibration, the meter relates this current to the concentration of DO.

CLARK ELECTRODE



MEMBRANE ELECTRODE METHOD

Calibration



Comparison with Winkler Titration

- 1. Fill two BOD bottles completely full of BOD dilution water, being very careful not to introduce air into either bottle.**
- 2. Analyze one bottle for D.O. using the Winkler titration.**
- 3. Insert the electrode into the second bottle, turn on the stirring mechanism, and wait for the reading to stabilize.**
- 4. Calibrate the meter to the D.O. value obtained in the titration.**
- 5. The meter is now ready for sample analysis**

MEMBRANE ELECTRODE METHOD

“AIR” Calibration

This procedure varies considerably among the various instrument models available. Therefore, the procedure must be obtained from the instrument manual, but the following points should be noted.



1. Where possible with the specific equipment being used, compensation should be made during calibration for both ambient temperature and local atmospheric pressure. This pressure should be determined using a reliable onsite barometer. The oxygen solubility table following this procedure may be used.
2. Carefully blot any water droplets from the membrane using a soft tissue.
3. During calibration, be sure the membrane is exposed to fresh air. Laying the electrode on the bench for calibration is usually adequate.
4. Complete the calibration as soon as possible before the electrode membrane begins to dry.
5. The temperature registered on the meter should be checked against a trusted thermometer often.
6. Daily calibration of the D.O. meter is required. Calibration should also be verified after every five or six sample measurements.
7. Assure sufficient sample flow across membrane surface during analysis to overcome erratic response.

	Onsite Barometric Pressure														
	Atm:	0.970	0.975	0.980	0.985	0.990	0.995	1.000	1.005	1.010	1.015	1.020	1.030	1.040	1.050
	mm Hg:	737	741	745	749	752	756	760	764	768	771	775	783	790	798
	Inch Hg:	29.02	29.17	29.32	29.47	29.62	29.77	29.92	30.07	30.22	30.37	30.52	30.82	31.12	31.42
Temperature° Celsius	15.00	9.78	9.83	9.88	9.93	9.98	10.03	10.08	10.13	10.19	10.24	10.29	10.39	10.49	10.60
	15.50	9.67	9.72	9.77	9.82	9.88	9.93	9.98	10.03	10.08	10.13	10.18	10.28	10.38	10.48
	16.00	9.57	9.62	9.67	9.72	9.77	9.82	9.87	9.92	9.97	10.02	10.07	10.17	10.27	10.37
	16.50	9.47	9.52	9.57	9.62	9.67	9.72	9.77	9.82	9.87	9.92	9.97	10.07	10.16	10.26
	17.00	9.37	9.42	9.47	9.52	9.57	9.62	9.66	9.71	9.76	9.81	9.86	9.96	10.06	10.16
	17.50	9.27	9.32	9.37	9.42	9.47	9.52	9.57	9.61	9.66	9.71	9.76	9.86	9.96	10.05
	18.00	9.18	9.23	9.27	9.32	9.37	9.42	9.47	9.51	9.56	9.61	9.66	9.76	9.85	9.95
	18.50	9.08	9.13	9.18	9.23	9.28	9.32	9.37	9.42	9.47	9.51	9.56	9.66	9.75	9.85
	19.00	8.99	9.04	9.09	9.13	9.18	9.23	9.28	9.32	9.37	9.42	9.47	9.56	9.65	9.75
	19.50	8.90	8.95	9.00	9.04	9.09	9.14	9.18	9.23	9.28	9.32	9.37	9.47	9.56	9.65
	20.00	8.81	8.86	8.91	8.95	9.00	9.05	9.09	9.14	9.18	9.23	9.28	9.37	9.45	9.56
	20.50	8.72	8.77	8.82	8.86	8.91	8.96	9.00	9.05	9.10	9.14	9.19	9.28	9.37	9.46
	21.00	8.64	8.68	8.73	8.78	8.82	8.87	8.91	8.96	9.01	9.05	9.10	9.19	9.28	9.37
	21.50	8.56	8.60	8.65	8.69	8.74	8.78	8.83	8.87	8.92	8.96	9.01	9.10	9.19	9.28
	22.00	8.47	8.52	8.56	8.61	8.65	8.70	8.74	8.79	8.83	8.88	8.92	9.01	9.10	9.19
	22.50	8.39	8.44	8.48	8.53	8.57	8.62	8.66	8.70	8.75	8.79	8.84	8.93	9.02	9.10
	23.00	8.31	8.36	8.40	8.44	8.49	8.53	8.58	8.62	8.67	8.71	8.75	8.84	8.93	9.02
	23.50	8.23	8.28	8.32	8.37	8.41	8.45	8.50	8.54	8.58	8.63	8.67	8.76	8.85	8.93
	24.00	8.16	8.20	8.24	8.29	8.33	8.37	8.42	8.46	8.50	8.55	8.59	8.68	8.76	8.85
	24.50	8.08	8.12	8.17	8.21	8.25	8.30	8.34	8.38	8.43	8.47	8.51	8.60	8.68	8.77
	25.00	8.01	8.05	8.09	8.13	8.18	8.22	8.26	8.31	8.35	8.39	8.43	8.52	8.60	8.69
	25.50	7.93	7.97	8.02	8.06	8.10	8.15	8.19	8.23	8.27	8.31	8.36	8.44	8.53	8.61
	26.00	7.86	7.90	7.94	7.99	8.03	8.07	8.11	8.15	8.20	8.24	8.28	8.36	8.45	8.53
	26.50	7.79	7.83	7.87	7.91	7.96	8.00	8.04	8.08	8.12	8.16	8.21	8.29	8.37	8.46
	27.00	7.72	7.76	7.80	7.84	7.88	7.93	7.97	8.01	8.05	8.09	8.13	8.22	8.30	8.38
	27.50	7.65	7.69	7.73	7.77	7.81	7.86	7.90	7.94	7.98	8.02	8.06	8.14	8.22	8.31
	28.00	7.58	7.62	7.66	7.70	7.75	7.79	7.83	7.87	7.91	7.95	7.99	8.07	8.15	8.23
	28.50	7.52	7.56	7.60	7.64	7.68	7.72	7.76	7.80	7.84	7.88	7.92	8.00	8.08	8.16
	29.00	7.45	7.49	7.53	7.57	7.61	7.65	7.69	7.73	7.77	7.81	7.85	7.93	8.01	8.09
	29.50	7.38	7.42	7.46	7.50	7.54	7.58	7.62	7.66	7.70	7.74	7.78	7.86	7.94	8.02
	30.00	7.32	7.36	7.40	7.44	7.48	7.52	7.56	7.60	7.64	7.68	7.72	7.79	7.87	7.95

NOTE: The first three lines are different units for the same pressure measurement.

CLARK ELECTRODE

Calibration

NOT
REQUIRED

NO



Temperature

Atmospheric Pressure


(Use a Reliable Barometer)



Luminescence D.O. Probe




Luminescence D.O. Probe




How Does LDO Work?

- A sensor is coated with a luminescent material.
- Blue light from an LED strikes the luminescent chemical on the sensor.
- The luminescent chemical instantly becomes excited.




How Does LDO Work?

- As the **excited chemical relaxes**, it **releases red light**.
- The red light is detected by a photo diode.
- The **time** it takes for the chemical to return to a relaxed state is measured




How Does LDO Work?

- When **oxygen** contacts the luminescent chemical, the intensity of the red light decreases
- The amount of time it takes for the material to relax is reduced




How Does LDO Work?

- The intensity of the red light is not what's being measured.
- What's being measured is the time it takes after excitation for red light to be given off.
 - Lifetime of luminescence



How Does LDO Work?

- A red LED is also present in the probe.
- Between flashes of the blue LED, a red LED of known intensity, is flashed on the sensor.
- The red LED acts as an **internal standard** (or reference) for a comparison to the red light given off by the luminescent chemical.



Why is this a Big Deal?

Reduced Maintenance

No membrane to replace

No more stretching of Teflon and worrying about air bubbles

No more punctured membranes

No electrolyte to foul or poison

No H₂S poisoning of the electrolyte

No anode or cathode

No cleaning of anodes

No more coating of electrodes

Why is this a Big Deal?

Reduced Maintenance

No membrane to replace

No more stretching of Teflon and worrying about air bubbles
No more punctured membranes

No electrolyte to foul or poison

No H₂S poisoning of the electrolyte

No anode or cathode

No cleaning of anodes
No more coating of electrodes

Why is this a Big Deal?

Reduced Maintenance

No membrane to replace

No more stretching of Teflon and worrying about air bubbles
No more punctured membranes

No electrolyte to foul or poison

No H₂S poisoning of the electrolyte

No anode or cathode

No cleaning of anodes
No more coating of electrodes

Why is this a Big Deal?

- **Frequent Calibration Not Required**
 - No anode to consume and no electrolyte to deplete means extremely stable measurements
 - Internal standard with Red LED
 - No interference from pH swings, wastewater chemicals, H₂S, or heavy metals

Why is this a Big Deal?

- **Accurate and Stable Readings**
 - With nothing to interfere with the readings, LDO produces more stable measurements for a longer time
- **Speed!**
 - Turn it on and it's running!
 - Response time of less than 30 seconds to 90%!

Why is this a Big Deal?

- **Simple Operation and Maintenance**
 - Only one replacement part
 - Inexpensive sensor cap is simple to replace quickly

NOW EPA APPROVED

Listed In
Federal Register
Vol. 77, No. 97
Friday, May 18, 2012

Listed as ASTM Method D888-09 (C)

Footnote 63 – Hach Method 10360
(BOD and cBOD)

Footnote 64 In-Situ Method 1002-8-2009
(Dissolved Oxygen Measurement by Optical Probe)

D.O. Meter



D.O. METER

Advantages

Saves Time

Continuous Monitoring

Less Chemical Interference

Portable

D.O. METER

Limitations

(Membrane Electrode)

Daily Calibration

Flow Past Membrane

Membrane May Foul

Requires Training

BOD₅

B.O.D.

Biochemical Oxygen Demand

**The Quantity of Oxygen Used in the
Biochemical Oxidation of Organic Material.**

Under:

Specified Time

5 Days


Specified Temperature

20⁰ C

Specified Conditions

**In the Dark
In the Presence
of Bacteria**

Measured Volume of Wastewater is Added to BOD Bottle.




Contains:

Organics

Ammonia

Bacteria

Dilution Water Added



Dilution Water Contains:

Nutrients

Oxygen


**Measure
D.O. Concentration**

Incubate 5 Days

Some Oxygen Used:


Respiration

Nitrification



**Measure
D.O. Concentration**

Measure Oxygen Loss (Demand)



D.O. Day 1 - D.O. Day 5

= Oxygen Demand

**(Of What Is In The
Bottle)**

BOD

REAGENTS

Distilled Water

High Quality

Free of Toxic Material

Free of Oxygen Demanding Substances

BOD

REAGENTS

Distilled Water

Phosphate Buffer

Magnesium Sulfate

Calcium Chloride

Ferric Chloride

Provide Essential Nutrients

Buffer pH



OTHER REAGENTS

DECHLORINATING AGENT

Sodium Sulfite - Na_2SO_3

NITRIFICATION INHIBITOR

CBOD

QUALITY CONTROL CHECK

Accuracy

Glucose - Glutamic Acid Solution

SAMPLE PRETREATMENT

Temperature

Near 20°C

pH

Between 6.5 and 7.5

(Adjust if > 8.5 or < 6.0 and seed)

Supersaturated D.O.

Agitate

Dechlorinate

Proper Amount of Sodium Sulfite

DECHLORINATION

100 mL of Sample

+ Potassium Iodide + Sulfuric Acid + Starch



**Titrate with Sodium
Sulfite**

**to Starch
Iodide
Endpoint**

DILUTION WATER

Distilled Water

plus

BUFFER

plus

NUTRIENTS

High Quality

No Toxics

No Organics

pH 7.2

**(Phosphorus
and Ammonia)**

Magnesium

Calcium

Iron


DILUTION WATER PREPARATION

NEEDED VOLUME for Each Day's Use

ADD NUTRIENTS (1 mL per Liter)

SATURATE WITH OXYGEN

**Shake
(small volume)
or
Draw Vacuum**



DILUTION WATER **PREPARATION**

NEEDED VOLUME

ADD NUTRIENTS

SATURATE WITH OXYGEN

STORE

In Incubator

ADD BUFFER

Day of Use

BOD PROCEDURE

DILUTE SAMPLE

Minimum Residual, 1.0 mg/L

Minimum Depletion, 2.0 mg/L

At Least Two Dilutions

Thoroughly Mix Sample

ADD NITRIFICATION INHIBITOR

If Required for CBOD



TCMP
0.10 gram/bottle
Two "shots"



BOD PROCEDURE

DILUTE SAMPLE

Minimum Residual, 1.0 mg/L

Minimum Depletion, 2.0 mg/L

At Least Two Dilutions

Thoroughly Mix Sample

ADD NITRIFICATION INHIBITOR

If Required for CBOD

Add SEED (Bacteria)

If Required

Disinfected Samples

Industrial Samples

Reference Samples

BOD PROCEDURE

Source of Seed (bacteria)

Settled Sewage

Primary Effluent



**Commercially
Available**



Adapted Seed

BOD PROCEDURE, (cont.)

FILL BOTTLE (with dilution water)

MEASURE INITIAL D.O.



STOPPER

(No Air Bubbles)

SEAL

(Water and Cover)



BOD PROCEDURE, (cont.)

INCUBATE

$20 \pm 1^{\circ}\text{C}$

5 Days \pm 6 hour



MEASURE FINAL D.O.

(wash bottles)

DILUTION WATER **BLANK**

CHECK ON QUALITY

MAXIMUM DEPLETION
0.2 mg/L

NOT USED IN CALCULATIONS

B.O.D. Calculation

**D.O. In - D.O. Out
= DEPLETION**

= Oxygen Demand of Diluted Sample

$$C_1 \times V_1 = C_2 \times V_2$$

$$\text{Sample BOD} \times \text{Sample Volume} = \text{BOD Diluted} \times \text{Diluted Volume}$$

$$\text{BOD Sample} = \frac{\text{BOD Diluted}}{\text{Sample Volume}} \times \text{Volume Diluted}$$

$$\text{BOD Sample} = \frac{\text{Depletion}}{\text{Sample Volume}} \times 300 \text{ mL}$$

B.O.D. mg/L =

$$\frac{\text{D.O. DEPLETION (mg/L)}}{\text{SAMPLE VOLUME (mL)}} \times 300 \text{ mL}$$

D.O. DEPLETION = D.O. Initial - D.O. 5-Day

Minimum Depletion - 2.0 mg/L

Minimum Residual - 1.0 mg/L

B.O.D. Example Problem

Calculate the B.O.D.:

Initial Sample D.O. = 8.1 mg/L
5-Day Sample D.O. = 2.1 mg/L
Vol. Of Sample in 300 mL Bottle = 60 mL

$$\text{B.O.D., mg/L} = \frac{\text{D.O. Depletion, mg/L}}{\text{Volume Sample, mL}} \times 300 \text{ mL}$$

$$= \frac{\text{Initial D.O. - Residual D.O.}}{\text{Volume Sample, mLs}} \times 300 \text{ mL}$$

$$= \frac{8.1 \text{ mg/L} - 2.1 \text{ mg/L}}{60 \text{ mL}} \times 300 \text{ mL}$$

$$= \frac{6.0 \text{ mg/L}}{60 \text{ mL}} \times 300 \text{ mL}$$

$$= 0.1 \text{ mg/L} \times 300 = 30 \text{ mg/L}$$

B.O.D. PRACTICE PROBLEM

Calculate the B.O.D. value to be reported for each of the samples below. Be sure to consider the minimum depletion and residual requirements.

<u>Sample 1</u>	<u>Dilution A</u>	<u>Dilution B</u>
mL Sample	15	30
Initial D.O., mg/L	8.1	8.2
5-Day D.O., mg/L	6.6	<u>4.2</u>
Depletion	1.5	4.0

$$\text{B.O.D., mg/L} = \frac{\text{D.O. Depletion, mg/L}}{\text{Volume Sample, mL}} \times 300 \text{ mL}$$

$$= \frac{4.0 \text{ mg/L}}{30 \text{ mL}} \times 300 \text{ mL}$$

$$= 40 \text{ mg/L}$$

B.O.D. PRACTICE PROBLEM

Calculate the B.O.D. value to be reported for each of the samples below. Be sure to consider the minimum depletion and residual requirements.

<u>Sample 2</u>	<u>Dilution A</u>	<u>Dilution B</u>
mL Sample	15	30
Initial D.O., mg/L	8.0	8.2
5-Day D.O., mg/L	<u>4.3</u>	<u>0.7</u>
Depletion	3.7	

$$\text{B.O.D., mg/L} = \frac{\text{D.O. Depletion, mg/L}}{\text{Volume Sample, mL}} \times 300 \text{ mL}$$

$$= \frac{3.7 \text{ mg/L}}{15 \text{ mL}} \times 300 \text{ mL}$$

$$= 74 \text{ mg/L}$$

B.O.D. PRACTICE PROBLEM

<u>Sample 3</u>	<u>Dilution A</u>	<u>Dilution B</u>
mL Sample	15	30
Initial D.O., mg/L	8.1	8.1
5-Day D.O., mg/L	<u>5.6</u>	<u>3.3</u>
Depletion	2.5	4.8

$$\text{BOD}_A = \frac{2.5 \text{ mg/L}}{15 \text{ mL}} \times 300 \text{ mL} = 50 \text{ mg/L}$$

$$\text{BOD}_B = \frac{4.8 \text{ mg/L}}{30 \text{ mL}} \times 300 \text{ mL} = 48 \text{ mg/L}$$

$$\frac{50 + 48}{2} = 49 \text{ mg/L}$$

BOD PROCEDURE

DILUTE SAMPLE

Minimum Residual, 1.0 mg/L

Minimum Depletion, 2.0 mg/L

At Least Two Dilutions

Thoroughly Mix Sample

ADD NITRIFICATION INHIBITOR

If Required for CBOD

Add SEED (Bacteria)

If Required

De-chlorinated Samples

Industrial Samples

Reference Samples



BOD PROCEDURE

Source of Seed (bacteria)

**Settled Sewage
Primary Effluent**





**Commercially
Available**



Adapted Seed

SEEDED BOD PROCEDURE





10 mL of seed caused 5.0 mg/L Depletion

1 mL would have caused $\frac{5.0 \text{ mg/L}}{10 \text{ mL}} = 0.5 \text{ mg/L}$

**So: 0.5 mg/L would be used by 1 mL of seed
in seeded sample**

SEEDED BOD

CALCULATION

$$\text{BOD} = \frac{\text{Depletion Caused by Sample}}{\text{Volume of Sample Used}} \times 300 \text{ mL}$$

$$\text{BOD}_{(\text{seeded})} = \frac{D_1 - D_2}{\text{Sample Volume}} \times 300 \text{ mL}$$

SEEDED BOD CALCULATION

$$\text{BOD}_{(\text{seeded})} = \frac{D_1 - D_2}{\text{Sample Volume}} \times 300 \text{ mL}$$

D_1 = Depletion Due To Sample and Seed
(Total Depletion in Bottle with Seed and Sample)

D_2 = D.O. Depletion Due to Just Seed
(Seed Depletion in Bottle with Seed and Sample)



$$\text{BOD}_{(\text{seeded})} = \frac{D_1 - D_2}{\text{Sample Volume}} \times 300 \text{ mL}$$

D_2 = D.O. Depletion Due to Just Seed
(Seed Depletion in Bottle with Seed and Sample)

D_2 Must Be Calculated
Based on Bottle with Just Seed

$$D_2 = \frac{\text{D.O. Depletion of a Sample of Seed}}{\text{Volume of Seed Used}}$$

In The Bottle with Just Seed

$$\text{BOD}_{(\text{seeded})} = \frac{D_1 - D_2}{\text{Sample Volume}} \times 300 \text{ mL}$$

D_2 = D.O. Depletion Due to Just Seed
(Seed Depletion in Bottle with Seed and Sample)

$$D_2 = \frac{\text{D.O. Depletion of a Sample of Seed}}{\text{Volume of Seed Used}}$$

In the Example

$$D_2 = \frac{5.0 \text{ mg/L}}{10 \text{ mL}} = 0.5 \text{ mg/L}$$

$$\text{BOD}_{(\text{seeded})} = \frac{D_1 - D_2}{\text{Sample Volume}} \times 300 \text{ mL}$$

$$\text{BOD}_{(\text{seeded})} =$$

$$\frac{3.7 \text{ mg/L} - 0.5 \text{ mg/L}}{100 \text{ mL}} \times 300 \text{ mL}$$


$$= \frac{3.2 \text{ mg/L}}{100 \text{ mL}} \times 300 \text{ mL}$$

$$= 9.6 \text{ mg/L}$$

D_2 = D.O. Depletion Due to Just Seed
(Seed Depletion in Bottle with Seed and Sample)

**What If More Than 1mL
Was Used for Seeding?**

$$D_2 = \frac{\text{D.O. Depletion of a Sample of Seed}}{\text{Volume of Seed Used}}$$


$$D_2 = \frac{5.0 \text{ mg/L}}{10 \text{ mL}} = 0.5 \text{ mg/L/mL}$$


Multiply by mL Used for Seeding

Example Seeded BOD Calculation.

Calculate the B.O.D. of a reference sample given that 1mL of seed material was used in the reference sample. (300 mL BOD bottles were used)

	<u>Seed Material</u>	<u>Reference Sample</u>
Volume Used	5 mL	5 mL
Initial D.O.	8.6 mg/L	8.7 mg/L
Final D.O.	5.8 mg/L	4.8 mg/L



Example Seeded BOD Calculation.

Calculate the B.O.D. of a reference sample given that 1mL of seed material was used in the reference sample. (300 mL BOD bottles were used)

	<u>Seed Material</u>	<u>Reference Sample</u>
Volume Used	5 mL	5 mL
Initial D.O.	8.6 mg/L	8.7 mg/L
Final D.O.	5.8 mg/L	4.8 mg/L

$$\text{BOD}_{\text{seeded}} = \frac{D_1 - D_2}{\text{Sample Volume, mL}} \times 300 \text{ mL}$$

D_1 = Depletion of Sample & Seed

$$D_1 = 8.7 \text{ mg/L} - 4.8 \text{ mg/L}$$

$$D_1 = 3.9 \text{ mg/L}$$

Example Seeded BOD Calculation.

Calculate the B.O.D. of a reference sample given that 1mL of seed material was used in the reference sample. (300 mL BOD bottles were used)

	<u>Seed Material</u>	<u>Reference Sample</u>
Volume Used	5 mL	5 mL
Initial D.O.	8.6 mg/L	8.7 mg/L
Final D.O.	5.8 mg/L	4.8 mg/L

$$\text{BOD}_{\text{seeded}} = \frac{D_1 - D_2}{\text{Sample Volume, mL}} \times 300 \text{ mL}$$

D_2 = Depletion of Seed
(In Seeded Sample)

$$D_2 = \frac{8.6 \text{ mg/L} - 5.8 \text{ mg/L}}{5 \text{ mL}}$$

$$D_2 = \frac{2.8 \text{ mg/L}}{5 \text{ mL}} = 0.56$$

Example Seeded BOD Calculation.

Calculate the B.O.D. of a reference sample given that 1mL of seed material was used in the reference sample. (300 mL BOD bottles were used)

	<u>Seed Material</u>	<u>Reference Sample</u>
Volume Used	5 mL	5 mL
Initial D.O.	8.6 mg/L	8.7 mg/L
Final D.O.	5.8 mg/L	4.8 mg/L

$$\text{BOD}_{\text{seeded}} = \frac{D_1 - D_2}{\text{Sample Volume, mL}} \times 300 \text{ mL}$$

$$\text{BOD} = \frac{3.9 \text{ mg/L} - 0.56 \text{ mg/L}}{5 \text{ mL}} \times 300 \text{ mL}$$

$$\text{BOD} = \frac{3.34 \text{ mg/L}}{5 \text{ mL}} \times 300 \text{ mL}$$

$$\text{BOD} = 200.4 \text{ mg/L}$$

Example #2 Seeded BOD Calculation.

Calculate the B.O.D. of an de-chlorinated effluent sample given that 1mL of seed material was used in the sample. (300 mL BOD bottles were used)

	<u>Seed Material</u>	<u>Sample</u>
Volume Used	10 mL	200 mL
Initial D.O.	8.3 mg/L	8.5 mg/L
Final D.O.	5.6 mg/L	3.7 mg/L

$$\text{BOD}_{\text{seeded}} = \frac{D_1 - D_2}{\text{Sample Volume, mL}} \times 300 \text{ mL}$$

D_1 = Depletion of Sample & Seed

$$D_1 = 8.5 \text{ mg/L} - 3.7 \text{ mg/L}$$

$$D_1 = 4.8 \text{ mg/L}$$

Example #2 Seeded BOD Calculation.

Calculate the B.O.D. of an de-chlorinated effluent sample given that 1mL of seed material was used in the sample. (300 mL BOD bottles were used)

	<u>Seed Material</u>	<u>Sample</u>
Volume Used	10 mL	200 mL
Initial D.O.	8.3 mg/L	8.5 mg/L
Final D.O.	5.6 mg/L	3.7 mg/L

$$\text{BOD}_{\text{seeded}} = \frac{D_1 - D_2}{\text{Sample Volume, mL}} \times 300 \text{ mL}$$

D_2 = Depletion of Seed
(In Seeded Sample)

$$D_2 = \frac{8.3 \text{ mg/L} - 5.6 \text{ mg/L}}{10 \text{ mL}}$$

$$D_2 = \frac{2.7 \text{ mg/L}}{10 \text{ mL}} = 0.27$$

Example #2 Seeded BOD Calculation.

Calculate the B.O.D. of an de-chlorinated effluent sample given that 1mL of seed material was used in the sample. (300 mL BOD bottles were used)

	<u>Seed Material</u>	<u>Sample</u>
Volume Used	10 mL	200 mL
Initial D.O.	8.3 mg/L	8.5 mg/L
Final D.O.	5.6 mg/L	3.7 mg/L

$$\text{BOD}_{\text{seeded}} = \frac{D_1 - D_2}{\text{Sample Volume, mL}} \times 300 \text{ mL}$$

$$\text{BOD} = \frac{4.8 \text{ mg/L} - 0.27 \text{ mg/L}}{200 \text{ mL}} \times 300 \text{ mL}$$

$$\text{BOD} = \frac{4.53 \text{ mg/L}}{200 \text{ mL}} \times 300 \text{ mL}$$

$$\text{BOD} = 6.8 \text{ mg/L}$$

BOD Practice Problems

1. Calculate the B.O.D. given:

D.O. IN = 7.0 mg/L

D.O. OUT (5-day) = 3.5 mg/L

Vol. Sample in B.O.D. bottle = 15 mL

2. Calculate the B.O.D. of a reference sample given that 1mL of seed material was used in the reference sample.
(300 mL BOD bottles were used)

	<u>Seed Material</u>	<u>Reference Sample</u>
Volume Used	9 mL	5 mL
Initial D.O.	8.7 mg/L	8.6 mg/L
Final D.O.	5.1 mg/L	5.0 mg/L

Work Calculations on Separate Paper
Answers Given on Next Slides

1. Calculate the B.O.D. given:

D.O. IN = 7.0 mg/L

D.O. OUT (5-day) = 3.5 mg/L

Vol. Sample in B.O.D. bottle = 15 mL

$$\text{BOD} = \frac{\text{Depletion, mg/L}}{\text{Sample Volume, mL}} \times 300 \text{ mL}$$

$$\text{BOD} = \frac{7.0 \text{ mg/L} - 3.5 \text{ mg/L}}{15 \text{ mL}} \times 300 \text{ mL}$$

$$\text{BOD} = \frac{3.5 \text{ mg/L}}{15 \text{ mL}} \times 300 \text{ mL}$$

$$\text{BOD} = 70 \text{ mg/L}$$

2. Calculate the B.O.D. of a reference sample given that 1mL of seed material was used in the reference sample.

(300 mL BOD bottles were used)

	<u>Seed Material</u>	<u>Reference Sample</u>
Volume Used	9 mL	5 mL
Initial D.O.	8.7 mg/L	8.6 mg/L
Final D.O.	5.1 mg/L	5.0 mg/L

$$\text{BOD}_{\text{seeded}} = \frac{D_1 - D_2}{\text{Sample Volume, mL}} \times 300 \text{ mL}$$

D_1 = Depletion of Sample & Seed

$$D_1 = 8.6 \text{ mg/L} - 5.0 \text{ mg/L}$$

$$D_1 = 3.6 \text{ mg/L}$$

D_2 = Depletion of Seed
(In Seeded Sample)

$$D_2 = \frac{8.7 \text{ mg/L} - 5.1 \text{ mg/L}}{9 \text{ mL}}$$

$$D_2 = \frac{3.6 \text{ mg/L}}{9 \text{ mL}} = 0.40$$

2. Calculate the B.O.D. of a reference sample given that 1mL of seed material was used in the reference sample.
(300 mL BOD bottles were used)

	<u>Seed Material</u>	<u>Reference Sample</u>
Volume Used	9 mL	5 mL
Initial D.O.	8.7 mg/L	8.6 mg/L
Final D.O.	5.1 mg/L	5.0 mg/L

$$\text{BOD}_{\text{seeded}} = \frac{D_1 - D_2}{\text{Sample Volume, mL}} \times 300 \text{ mL}$$

D_1 = Depletion of Sample & Seed

D_2 = Depletion of Seed

$$\text{BOD} = \frac{3.6 \text{ mg/L} - 0.40 \text{ mg/L}}{5 \text{ mL}} \times 300 \text{ mL}$$

$$\text{BOD} = \frac{3.2 \text{ mg/L}}{5 \text{ mL}} \times 300 \text{ mL}$$

$$\text{BOD} = 192 \text{ mg/L}$$

Dissolved Oxygen and Biochemical Oxygen Demand Analyses



**Prepared By
Michigan Department of Environmental Quality
Operator Training and Certification Unit**